

## LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

Ref. [MEEMER03\\_O](#)

**EMERGENCY'S™**



13/01/2026 23:20:48 – For information purposes only, can be subject to change

## LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

Ref. MEEMER03\_O

**EMERGENCY'S™**

**Intended to transport the most fragile + relief materials**

### Description

- Four compartments for interior supports and high-capacity compartment.
- Transparent pouch in the inner flap, designed to see and organize the material easily, thanks to the rubber-elastic that guarantees a correct fixation-
- Five outer pockets, two side pockets, one on the back and one on the top and a front pocket that can be divided internally thanks to rubber pockets and bands
- Base of the waterproof bag
- High resistance zips
- Reflective strips-
- Reinforced handles and shoulder straps
- External card holder
- Dimensions: 55 x 31.5 x 33 cm, capacity 38 L
- Internal dimensions: 40 x 20 x 24 cm
- Material 600D Polyester
- Weight 1.76 kg
- Sold empty
- Available in : Red (MEEMER03\_R), Blue (MEEMER03\_B) and Orange (MEEMER03\_O)

## LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

**EMERGENCY'S™**

Ref. [MEEMER03\\_O](#)

### Features

|                             |               |                         |           |
|-----------------------------|---------------|-------------------------|-----------|
| CE Marking                  | NONE          | Medical device          | No        |
| Bar code                    | 3661040202891 | HS Code                 | 420299000 |
| Warranty (months)           | 24            | Country                 | China     |
| DIMENSIONS AND WEIGHT       |               |                         |           |
| Packed unit / Width (cm)    | 28            | Bare Unit / Width (cm)  | 31.5      |
| Packed unit / Height (cm)   | 20            | Bare Unit / Height (cm) | 33        |
| Packed unit / Length (cm)   | 50            | Bare Unit / Length (cm) | 55        |
| Packed unit / Weight        | 1.8           | Bare Unit / Weight (kg) | 1.76      |
| SPECIFIC STORAGE CONDITIONS |               |                         |           |
| Minimum temperature         | NA            | Max temperature         | NA        |
| Minimum humidity            | NA            | Maximum humidity        | NA        |
| OTHER FEATURES              |               |                         |           |
| Colors                      | Orange        | Packed unit / Quantity  | 1         |