

PULSE STETHOSCOPE / RED

Ref. 507522




CE



14/01/2026 21:36:38 – For information purposes only, can be subject to change

PULSE STETHOSCOPE / RED

Ref. 507522



Description

For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaura with built-in spring.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009281	HS Code	90181990
Warranty (months)	36		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.22

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Navy blue	Packed unit / Quantity	1
--------	-----------	------------------------	---

14/01/2026 21:36:38 – For information purposes only, can be subject to change