



Ref. MONI011E



CE



14/01/2026 21:35:07 – For information purposes only, can be subject to change

=====

Ref. MONI011E

=====

Description

=====

Features

CE Marking	CE	Medical device	Class I
Bar code	6948309002582	HS Code	9018193090
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Packed unit / Height (cm)	2
Packed unit / Length (cm)	18	Packed unit / Weight	0.15
Bare Unit / Weight (kg)	0.1		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES