

## REUSABLE NIBP FOR VITALIMATE MINI

Ref. MONI011



# REUSABLE NIBP FOR VITALIMATE MINI

Ref. MONI011

Vitalimate MINI reusable NIBP cuffs

## Description

The Vitalimate MINI NIBP cuffs are comfortable and come in all sizes for all types of patients.

- Reusable adult NIBP cuff (27-35cm) : ref MONI011
- Reusable large adult NIBP cuff (34-43cm) : ref MONI011A
- Small reusable adult NIBP cuff (20.5-28 cm): ref MONI011G
- Reusable pediatric NIBP cuff (14-21.5 cm): ref MONI011C
- Reusable neonatal NIBP cuff (10-15 cm): ref MONI011D
- Reusable adult NIBP thigh cuff (42-54 cm): ref MONI011B

Disposable neonatal cuffs are also available for infants.

## Features

| CE Marking        | CE            | Medical device | Class I    |
|-------------------|---------------|----------------|------------|
| Bar code          | 6948309002421 | HS Code        | 9018193090 |
| Warranty (months) | 6             | Country        | China      |

### DIMENSIONS AND WEIGHT

|                           |       |                           |      |
|---------------------------|-------|---------------------------|------|
| Packed unit / Width (cm)  | nc    | Packed unit / Height (cm) | nc   |
| Packed unit / Length (cm) | nc    | Packed unit / Weight      | 0.08 |
| Bare Unit / Weight (kg)   | 0.075 |                           |      |

### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -20 | Max temperature  | 55 |
| Minimum humidity    | 10  | Maximum humidity | 93 |

### OTHER FEATURES