

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520620



CE
0459



17/01/2026 17:51:56 – For information purposes only, can be subject to change

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520620



=====

Description

HIGH ACCURACY

- Single tubing with connectic with screw
- "Watchmaking" mechanism

SHOCK RESISTANT

- Non deformable metal casing
- Shock proof foam manometer protector

COMFORT

- Ambidextrous thanks to the central decompression valve
- Large pump for quick inflation
- SPENGLER Velcro fastening

CUFF

- Coton or or nylon cuff with velcro fastening
- Sizes: **XXS** (10-15 cm), **XS** (15-22 cm), **S** (20-28 cm), **M** (26-34 cm), **L** (32-42 cm) ou **XL** (42-52 cm)
- Available in blue, grey and Cartoon
- Available in Clinic® to prevent nosocomial infection
- Comes with a transport pouch

Adult (M), cotton, velcro, navy: ref 520627

Adult (M), cotton, velcro, grey: ref 520621

Adult (L), nylon, velcro, grey: ref 520620

Child/Small adult (S), cotton, velcro, motifs: ref 520622

Multicuff (L,M,S), velco: ref 520630

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520620



Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446024796	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.736
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		