

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520626



CE  
0459



## LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520626



=====

### Description

#### HIGH ACCURACY

- Single tubing with connectic with screw
- "Watchmaking" mechanism

#### SHOCK RESISTANT

- Non deformable metal casing
- Shock proof foam manometer protector

#### COMFORT

- Ambidextrous thanks to the central decompression valve
- Large pump for quick inflation
- SPENGLER Velcro fastening

#### CUFF

- Coton or or nylon cuff with velcro fastening
- Sizes: **XXS** (10-15 cm), **XS** (15-22 cm), **S** (20-28 cm), **M** (26-34 cm), **L** (32-42 cm) ou **XL** (42-52 cm)
- Available in blue, grey and Cartoon
- Available in Clinic® to prevent nosocomial infection
- Comes with a transport pouch

Adult (M), cotton, velcro, navy: ref 520627

Adult (M), cotton, velcro, grey: ref 520621

Adult (L), nylon, velcro, grey: ref 520620

Child/Small adult (S), cotton, velcro, motifs: ref 520622

Multicuff (L,M,S), velco: ref 520630

## LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520626



### Features

|                             |               |                           |                               |
|-----------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                  | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                    | 3700446042479 | HS Code                   | 90189010                      |
| Warranty (months)           | 36            | Country                   | France                        |
| DIMENSIONS AND WEIGHT       |               |                           |                               |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc                            |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.617                         |
| SPECIFIC STORAGE CONDITIONS |               |                           |                               |
| Minimum temperature         | -20           | Max temperature           | 70                            |
| Minimum humidity            | 85            | Maximum humidity          | 85                            |
| OTHER FEATURES              |               |                           |                               |
| Packed unit / Quantity      | 1             |                           |                               |