

## VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



**Ref. 330012**

### Description

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

Adult (M) cotton, velcro, grey: 330011

Adult (M) cotton, velcro, red : 330013

Adult (M) navy, with straps: 330112

Adult (M) red, with straps: 330113

Adult (L) cotton, velcro, navy: 330002

Child/small adult (S), cotton, velcro, grey: 330021

### Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | NONE          | Medical device | No         |
| Bar code          | 3700446028541 | HS Code        | 9018901000 |
| Warranty (months) | 6             | Country        | France     |

#### DIMENSIONS AND WEIGHT

|                           |      |                           |      |
|---------------------------|------|---------------------------|------|
| Packed unit / Width (cm)  | nc   | Packed unit / Height (cm) | nc   |
| Packed unit / Length (cm) | nc   | Packed unit / Weight      | 0.19 |
| Bare Unit / Weight (kg)   | 0.19 |                           |      |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|