

Ref. 336403



## Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

19/01/2026 21:48:24 – For information purposes only, can be subject to change

Ref. 336403



## Features

|                             |               |                           |            |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking                  | CE            | Medical device            | Class I    |
| Bar code                    | 3700446008819 | HS Code                   | 9018901000 |
| Warranty (months)           | 6             | Country                   | Indonesia  |
| DIMENSIONS AND WEIGHT       |               |                           |            |
| Packed unit / Width (cm)    | 18            | Packed unit / Height (cm) | 4          |
| Packed unit / Length (cm)   | 21            | Packed unit / Weight      | 0.272      |
| SPECIFIC STORAGE CONDITIONS |               |                           |            |
| Minimum temperature         | NA            | Max temperature           | NA         |
| Minimum humidity            | NA            | Maximum humidity          | NA         |
| OTHER FEATURES              |               |                           |            |
| Packed unit / Quantity      | 1             |                           |            |