

Ref. 513250



CE
0459



Description

Efficiency and high-performance improved with the Clinic® version. This version of the Lian® Nano is lightweight and easy to handle and efficient thanks to its Clinic® cuff, suitable for liquid decontamination. Ideal for septic isolation or for preventing cross-contamination. Supplied in a nylon pouch.

19/01/2026 21:22:59 – For information purposes only, can be subject to change

Ref. 513250



Features

CE Marking	CE0459	Medical device	No
Bar code	3700446024710	HS Code	90189010
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.565

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---