

Ref. 513252



CE  
0459



## Description

Efficiency and high-performance improved with the Clinic® version. This version of the Lian® Nano is lightweight and easy to handle and efficient thanks to its Clinic® cuff, suitable for liquid decontamination. Ideal for septic isolation or for preventing cross-contamination. Supplied in a nylon pouch.

19/01/2026 21:17:06 – For information purposes only, can be subject to change

Ref. 513252



## Features

CE Marking	CE0459	Medical device	No
Bar code	3700446024727	HS Code	90189010
Warranty (months)	24	Country	France

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.503

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---