


AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS106

DDM




CE



20/01/2026 22:55:08 – For information purposes only, can be subject to change

AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS106



Those lower limb cuffs are autoclavable!

Description

Lower limb cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

Dimensions of the bladder in cm

| | |
|--------------------------|--------------|
| Lower limb cuff infant | 35.5 x 6 cm |
| Lower limb cuff child | 53 x 6 cm |
| Lower limb cuff adult | 69 x 9 cm |
| Lower limb cuff adult L | 81 x 9 cm |
| Lower limb cuff adult XL | 85.5 x 12 cm |

AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. [GCS106](#)

Features

| | | | |
|-------------------|---------------|----------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3700468203681 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | France |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|-----|------------------|----|
| Minimum temperature | -10 | Max temperature | 40 |
| Minimum humidity | 30 | Maximum humidity | 40 |

OTHER FEATURES