

VITALIMATE MINI PATIENT MONITOR PNI+SPO2 SPENGLER



Ref. MON002



CE
0123



22/01/2026 20:34:56 – For information purposes only, can be subject to change

VITALIMATE MINI PATIENT MONITOR PNI+SPO2 SPENGLER



Ref. MON002

Vitalimate MINI the polyvalent portable monitor for all health care settings

Description

Vitalimate Mini measures SpO2, PR and NIBP. It is adapted for adult, pediatric and neonat patients use.

- High definition LCD color display 3,5" (480x320 pixels). Can display figures and plethysmographic waveform
- Un socle double fonction. Il permet une recharge du moniteur pour 12 heures d'utilisation continue ou 24 heures si laissé sur le socle. C'est aussi une base de rangement, gardant le moniteur prêt à l'emploi.
- Alarmes sonores et visuelles paramétrables.
- Format compacte, facile à transporter: 320g, dimensions 175 x 85 x 32 mm

Contenu de la boîte: un moniteur, le socle, un capteur SpO2 rigide adulte, un brassards PNI adulte, un chargeur USB et un manuel.

D'autres brassards PNI sont disponibles à l'achat séparément:

- Brassard PNI réutilisable adulte (27-35cm) : ref MONI011
- Brassard PNI réutilisable adulte large (34-43cm) : ref MONI011A
- Brassard PNI réutilisable petit adulte large (20.5-28cm) : ref MONI011G
- Cuissard PNI réutilisable adulte (42-54cm) : ref MONI011B
- Brassard PNI réutilisable enfant (18-26cm) : ref MONI011C
- Brassard PNI réutilisable néonatal (10-19cm) : ref MONI011D
- Lot de 3 brassards néonatales jetables (3 tailles): ref MONI011E

VITALIMATE MINI PATIENT MONITOR PNI+SPO2 SPENGLER



Ref. MON002

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	6948309002360	HS Code	9018193010
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	9	Packed unit / Height (cm)	33
Packed unit / Length (cm)	22	Packed unit / Weight	1.43
Bare Unit / Weight (kg)	0.66		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES