



Ref. TI50003S



CE
0459



📅 22/01/2026 21:42:26 – For information purposes only, can be subject to change

=====

Ref. TI50003S



=====

Description

=====

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215516	HS Code	90189010
Warranty (months)	36	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	20
Packed unit / Length (cm)	13	Packed unit / Weight	0.485
Bare Unit / Weight (kg)	0.325		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES