


Ref. TS53000S



CE
0459



📅 22/01/2026 21:38:58 – For information purposes only, can be subject to change

=====

Ref. TS53000S



=====

Description

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215486	HS Code	90189010
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	25
Packed unit / Length (cm)	16	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

 22/01/2026 21:38:58 – For information purposes only, can be subject to change