



=====

**Ref. GCCMS105**

=====

## Description

=====

## Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | NONE          | Medical device | Class I  |
| Bar code          | 3700468208457 | HS Code        | 90189084 |
| Warranty (months) | 24            | Country        | France   |

### DIMENSIONS AND WEIGHT

|                           |    |                           |     |
|---------------------------|----|---------------------------|-----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc  |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.1 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES