



=====
Ref. GCCMS106
=====

Description
=====

Features

| | | | |
|-------------------|---------------|----------------|----------|
| CE Marking | NONE | Medical device | Class I |
| Bar code | 3700468208464 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | France |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

27/01/2026 23:52:13 – For information purposes only, can be subject to change