


BAG EIFFEL, DE BOISSY

Ref. MDEIFF01

De Boissy
— DEPUIS 1982



📅 27/01/2026 23:53:01 – For information purposes only, can be subject to change

BAG EIFFEL, DE BOISSY

Ref. MDEIFF01

De Boissy
— DEPUIS 1982**The chic and simple large volume bag, easy to organize****Description**

Chic and simple at the same time, the large EIFFEL bag recalls the shape of the traditional leather bag of doctors. Thanks to its very large bidirectional opening and multiple pockets/compartments, it makes it easy to store and reach all the care equipment and diagnostic instruments that health professionals need. Of a mottled gray, the polyester fabric is resistant, durable, water-repellent and therefore washable, it is accompanied by a reflective strip for safety, and YKK® slides for their quality.

- The large volume interior main space includes 8 pockets, 2 removable velcro dividers, 1 zipped mesh pouch. With a very large bidirectional zip opening, it gives very easy access to diagnostic instruments and care equipment
- A ventral pocket (40x20 cm) opens completely with on one side 2 zipped storage nets and 2 flat pockets, and on the other side several elastic loops that allow to store small care supplies: pliers, scissors, pen lamp, thermometer, bandages, pens, etc...
- Four external zipped pockets, one for A4 documents and another on the side with elastic to slide the bottle of hydroalcoholic gel which remains easily accessible at any time
- Padded and solid handles, a padded and adjustable shoulder strap - A name tag to identify the bag
- A non-slip bottom to reduce friction and wear
- Dimensions w x h x d: 40 x 35 x 24 cm- Weight: 1.3 kg
- Maximum recommended content: 6 kg / 24 liters
- Sold empty

BAG EIFFEL, DE BOISSY

Ref. MDEIFF01

De Boissy
 — DEPUIS 1982

Features

| | | | |
|-------------------|---------------|----------------|-----------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3661040218939 | HS Code | 420229000 |
| Warranty (months) | 12 | Country | China |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|----|-------------------------|------|
| Packed unit / Width (cm) | 37 | Bare Unit / Width (cm) | 35 |
| Packed unit / Height (cm) | 10 | Bare Unit / Height (cm) | 24 |
| Packed unit / Length (cm) | 45 | Bare Unit / Length (cm) | 40 |
| Packed unit / Weight | nc | Bare Unit / Weight (kg) | 1.29 |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|