

**ELITE BAGS-EXTREME-  
POLYESTER/WATERPROOF-RED**

Ref. MEEXTR02\_R

**ELITE BAGS**



GARANTIE  
**2**  
ANS  
2 YEARS  
GUARANTEE



27/01/2026 23:48:51 – For information purposes only, can be subject to change

## ELITE BAGS-EXTREME- POLYESTER/WATERPROOF-RED



Ref. MEEXTR02\_R

### Description

For users of emergency services specialized in first aid.

Interior totally adaptable to the needs of the user thanks to extractable velcro separators.

Elastic bands in all interior flat surfaces.

Reflective strips night and day. 2 external pockets of medium size, 2 external pockets of large capacity, with elastic bands.

4 extractable compartments with transparent frontal part differentiated according to the standardized color code (blue: respiratory system, red: circulatory system, green: various surgical equipment, yellow: pediatric equipment).

Isothermal amp maker with 27 elastics.

Shoulder strap, handles that surround the bag for better resistance, concealed backpack straps and banana.

Material: 1000D polyamide fabric

Colour: Red.

Weight: 2.86 kg

Outdoor measures: 47 x 22 x 28 cm / 23.23 l. Maximum recommended capacity: 12 Kg

# ELITE BAGS-EXTREME-POLYESTER/WATERPROOF-RED



Ref. [MEEXTR02\\_R](#)

## Features

CE Marking	NONE	Medical device	No
Bar code	8436564540300	HS Code	4202129900
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	23	Bare Unit / Width (cm)	22
Packed unit / Height (cm)	39	Bare Unit / Height (cm)	28
Packed unit / Length (cm)	43	Bare Unit / Length (cm)	47
Packed unit / Weight	3.8	Bare Unit / Weight (kg)	2.86
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 27/01/2026 23:48:51 – For information purposes only, can be subject to change