

COMBAT APPLICATION TOURNIQUET, BLACK COLOR

Ref. STGARR01



📅 27/01/2026 23:51:37 – For information purposes only, can be subject to change

COMBAT APPLICATION TOURNIQUET, BLACK COLOR

Ref. STGARR01



Tourniquet very used by the military and rescuers.

Description

-

Strong and resistant •

Ready to use to intervene quickly in case of hemorrhage.

- Easy to apply, it can be used with one hand, and it allows rescuers to quickly free themselves to go see other victims.

-

Suitable for arms and leg (including thigh)

-

Dimensions: 90 x 3.8 cm•

A label is present to write the time of the withers inserted.

- Black color, it also exists in Orange

**COMBAT APPLICATION TOURNIQUET,
BLACK COLOR****Ref. STGARR01****Features**

CE Marking	CE	Medical device	Class I
Bar code	3361040000205	HS Code	300610900
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.078
Bare Unit / Weight (kg)	0.078		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---