

TEMPO EAR PRO + - EAR THERMOMETER (RECHARGEABLE)

Ref. THE003



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



27/01/2026 23:14:00 – For information purposes only, can be subject to change

TEMPO EAR PRO + - EAR THERMOMETER (RECHARGEABLE)

Ref. THE003



Performance, hygiene and comfort for healthcare professionals

Description

Spengler's Tempo Ear Pro is the infrared ear thermometer designed for demanding medical environments. Fast, reliable and hygienic, it makes daily temperature taking easier, whether in the office, hospital or clinic.

- Ultra-fast measurement: less than 3 seconds
- Probe preheated to 36°C: patient comfort
- Infrared technology: reliable, even when moving
- Probe cover ejection: maximum hygiene without contact
- Memory of last 30 time-stamped measurements
- Automatic sleep mode after 60 seconds
- Pulse counter
- Lanyard
- Rechargeable

Included with Tempo Ear Pro + :

- 1 Tempo Ear Pro + ear thermometer
- 1 storage and charging cradle
- 20 probe covers
- 1 attachment cord
- 1 user manual
- Thermometer calibration and verification certificate with batch number

Also available : Tempo EAR PRO with batteries instead of battery.

TEMPO EAR PRO + - EAR THERMOMETER (RECHARGEABLE)



Ref. THE003

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6970392219448	HS Code	9025190090
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	5	Packed unit / Height (cm)	3
Packed unit / Length (cm)	14	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.125		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	50
Minimum humidity	15	Maximum humidity	93

OTHER FEATURES

 27/01/2026 23:14:00 – For information purposes only, can be subject to change