

**Ref. 334150**



CE



## Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

 28/01/2026 22:12:45 – For information purposes only, can be subject to change

=====

**Ref. 334150**



## Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE            | Medical device | Class I    |
| Bar code          | 3700446026622 | HS Code        | 9018901000 |
| Warranty (months) | 6             | Country        | Indonesia  |

### DIMENSIONS AND WEIGHT

|                           |    |                           |       |
|---------------------------|----|---------------------------|-------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.335 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|