


LUXASCOPE AURIS LED 2.5 V - LUXAMED

Ref. A1.416.314

LUXAMED®



CE



LUXASCOPE AURIS LED 2.5 V - LUXAMED

Ref. A1.416.314

LUXAMED®

Innovative otoscope in pocket size

Description

- Illumination: LED-Ring with 6 LEDs (concentric),
no fibre optic, field of vision is not restricted
- Maintenance-free, because no lamp replacement necessary
 - Can be dimmed
 - Automatic switchoff after 3 minutes
 - Life expectancy of LEDs: approx. 100.000 hours
by constant current drive
 - Intensity of illumination: approx. 10.000 Lux
 - Colour temperature: approx. 4.200 K
 - Optimized design of head
 - Pivoting 3-fold magnifying glass with locking function
 - Insufflation port for pneumatic test
 - Two-component handle (aluminium / fibre-glass reinforced),
also combinable with LuxaScope LED Dermatoscope- and
Ophthalmoscope head
 - Electronic ON/OFF push-button switch
 - Power supply: 2x AAA Alkaline batteries, inclusive
 - Incl. 20 disposable ear funnels, grey,
each 10 x 2.5 mm and 10 x 4.0 mm

LUXASCOPE AURIS LED 2.5 V - LUXAMED

Ref. A1.416.314

LUXAMED®

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960034	HS Code	9018908400
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	11	Packed unit / Height (cm)	4
Packed unit / Length (cm)	24	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.12		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Gray
--------	------