

**Ref. BT101S**



CE



28/01/2026 21:15:01 – For information purposes only, can be subject to change

=====

**Ref. BT101S**

=====

## Description

=====

## Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE            | Medical device | Class I    |
| Bar code          | 3700468215660 | HS Code        | 9018901000 |
| Warranty (months) | 12            | Country        | France     |

### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |     |
|---------------------|----|------------------|-----|
| Minimum temperature | 10 | Max temperature  | 70  |
| Minimum humidity    | 10 | Maximum humidity | 100 |

### OTHER FEATURES