

Ref. BT103S



28/01/2026 21:59:26 – For information purposes only, can be subject to change

=====

Ref. BT103S

=====

Description

=====

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215684	HS Code	9018901000
Warranty (months)	12	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES