

TRABECULECTOMY SET

Ref. **IBOI418**



Description

This box contains:

- *Vannas corneal scissors, 10 cm, straight, soft
- Barraquer cannula, curved, olive-shaped, for cyclodialysis
- Scalpel handle no. 3 -IS
- Stainless steel box, 20 x 10 x 5 cm -IS
- Micro Halstead forceps, 10 cm, straight, A/G, Holtex
- *Barraquer Fixing Forceps, 1 x 2 oblique claws
- Elschmig Fixing Forceps, A/G, oblique
- *Bonn Corneal Forceps
- Paufigue Suture Forceps
- Barraquer Needle Holder, curved, 13 cm, universal jaws
- *Castroviejo Needle Holder, foam, 14 cm, curved
- Jaeger Spatula

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040082813	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.743
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		