

DIAGNOSTIC PENLIGHT LED - LUXAMED

Ref. D1.211.112

LUXAMED®



05/02/2026 00:06:16 – For information purposes only, can be subject to change

DIAGNOSTIC PENLIGHT LED - LUXAMED

Ref. D1.211.112

LUXAMED®

Diagnostic Penlight LED OF Luxamed is resistant and lasts over time thanks to LED technology!

This technology provides a shadow-free field of vision and optimum illumination in the area of application. The LEDs are powered by a constant current, guaranteeing considerably longer life.

Description

Penlight for optimal eye and throat examinations, with LED light approved for eye examinations.

- Top quality ABS plastic housing
- Weight 31.5 g
- Life time LED 100.000 h
- Luminous intensity 18.000 – 22.000 mcd
- Colour temperature: approx. 7000 Kelvin
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

White (Ref : D1.211.912)

Blue (Ref : D1.211.212)

Grey (Ref : D1.211.312)

Yellow (Ref : D1.211.512)

Kids (Ref : D1.311.912)

Black (Ref : D1.211.112)

Red (Ref : D1.211.412)

Green (Ref : D1.211.612)

DIAGNOSTIC PENLIGHT LED - LUXAMEDRef. D1.211.112**LUXAMED®****Features**

CE Marking	CE	Medical device	Class I
Bar code	04260381961673	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.35
Bare Unit / Weight (kg)	0.04		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Black
--------	-------