

LIAN METAL / GREY NYLON CUFF / OBESE (XL)

Ref. 521310



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE


GARANTIE
3
3 ANS
3 YEARS
GUARANTEE

📅 05/02/2026 23:52:23 – For information purposes only, can be subject to change

LIAN METAL / GREY NYLON CUFF / OBESE (XL)



Ref. 521310



Description

Large dial, ambidextrous model (simply rotate to change hands). Metal casing with large dial for easy reading from all angles. Soft PVC manometer cover. Ergonomic spoon handle and inflation bulb. Cuffs with screw-on metal quick coupling. Supplied in a nylon pouch.

Available in :

- Adult (M), cotton, velcro, red (521317)
- Adult (M), cotton, velcro, grey (521311)
- Adult large (L), cotton, velcro, grey (521310)
- Child/Small Adult (S), nylon, velcro, grey (521316)
- Multibrassard (S, M, L) velcro (521330)

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446024864	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.782
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 05/02/2026 23:52:23 – For information purposes only, can be subject to change