

TEMPO WATERPROOF DIGITAL THERMOMETER SPENGLER

Ref. THE001



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

05/02/2026 23:10:42 – For information purposes only, can be subject to change

TEMPO WATERPROOF DIGITAL THERMOMETER SPENGLER

Ref. THE001

TEMPO is a simple and precise digital thermometer

Description

TEMPO is a waterproof digital thermometer allowing 3 different uses for taking temperature: the oral, the rectal and the armpit one.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6970392217642	HS Code	9025199090
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.01
Bare Unit / Weight (kg)	0.01		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES