

# STERILIZING STAINLESS STEEL BOX 25X12X6 CM

Ref. AD25126



09/02/2026 13:17:43 – For information purposes only, can be subject to change

# STERILIZING STAINLESS STEEL BOX 25X12X6 CM

Ref. [AD25126](#)



**Convenient for storing surgical instruments**

## Description

- Sold with lid to protect instruments
- 
- Decontamination with detergent and disinfectant -
- 
- Sterilizable by autoclave (134°-18 minutes)
- 
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- 
- Available in different sizes

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040013770	HS Code	90189084
Warranty (months)	12	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Bare Unit / Width (cm)	12
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	27	Bare Unit / Length (cm)	25
Packed unit / Weight	0.744	Bare Unit / Weight (kg)	0.726

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 09/02/2026 13:17:43 – For information purposes only, can be subject to change