

## S/S BOX 18 X 8 X 4 CM

Ref. AI18084



09/02/2026 13:26:34 – For information purposes only, can be subject to change

## S/S BOX 18 X 8 X 4 CM

Ref. [AI18084](#)



**Convenient for storing surgical instruments**

### Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)- In the " ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately
- Available in different sizes

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019482	HS Code	90189084
Warranty (months)	12	Country	Pakistan

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	19	Bare Unit / Width (cm)	8
Packed unit / Height (cm)	10	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	18
Packed unit / Weight	0.265	Bare Unit / Weight (kg)	0.256

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 09/02/2026 13:26:34 – For information purposes only, can be subject to change