

S/S BOX 25 X 10 X 5 CM

Ref. AI25105



09/02/2026 13:23:00 – For information purposes only, can be subject to change

S/S BOX 25 X 10 X 5 CM

Ref. [AI25105](#)



Ideal for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

Features

| | | | |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040019543 | HS Code | 90189084 |
| Warranty (months) | 12 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 10 | Bare Unit / Width (cm) | 10 |
| Packed unit / Height (cm) | 5 | Bare Unit / Height (cm) | 5 |
| Packed unit / Length (cm) | 25 | Bare Unit / Length (cm) | 25 |
| Packed unit / Weight | 0.563 | Bare Unit / Weight (kg) | 0.563 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

📅 09/02/2026 13:23:00 – For information purposes only, can be subject to change