

**S/S BOX 40 X 20 X 9 CM**

**Ref. AI42188**



09/02/2026 12:58:19 – For information purposes only, can be subject to change

## S/S BOX 40 X 20 X 9 CM

Ref. [AI42188](#)



**Convenient for storing surgical instruments**

### Description

- Sold with lid to protect instruments
- 
- Decontamination with detergent and disinfectant -
- 
- Sterilizable by autoclave (134°-18 minutes)
- 
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- 
- Available in different sizes

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019758	HS Code	90189084
Warranty (months)	12	Country	Pakistan

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	20	Bare Unit / Width (cm)	20
Packed unit / Height (cm)	9	Bare Unit / Height (cm)	9
Packed unit / Length (cm)	40	Bare Unit / Length (cm)	40
Packed unit / Weight	1.75	Bare Unit / Weight (kg)	1.749

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

 09/02/2026 12:58:19 – For information purposes only, can be subject to change