

**S/S BOX 50 X 20 X 12 CM**

**Ref. AI50201**



09/02/2026 13:17:04 – For information purposes only, can be subject to change

# S/S BOX 50 X 20 X 12 CM

Ref. [AI50201](#)



**Convenient for storing surgical instruments**

## Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019772	HS Code	90189084
Warranty (months)	12	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	50	Bare Unit / Width (cm)	20
Packed unit / Height (cm)	12	Bare Unit / Height (cm)	12
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	50
Packed unit / Weight	2.263	Bare Unit / Weight (kg)	2.263

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 09/02/2026 13:17:04 – For information purposes only, can be subject to change