


# AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS106

**DDM** 



CE



09/02/2026 12:52:34 – For information purposes only, can be subject to change

# AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS106



Those lower limb cuffs are autoclavable!

## Description

Lower limb cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

### Dimensions of the bladder in cm

|                          |              |
|--------------------------|--------------|
| Lower limb cuff infant   | 35.5 x 6 cm  |
| Lower limb cuff child    | 53 x 6 cm    |
| Lower limb cuff adult    | 69 x 9 cm    |
| Lower limb cuff adult L  | 81 x 9 cm    |
| Lower limb cuff adult XL | 85.5 x 12 cm |

# AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. GCS106

## Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3700468203681 | HS Code        | 90189084 |
| Warranty (months) | 24            | Country        | France   |

### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -10 | Max temperature  | 40 |
| Minimum humidity    | 30  | Maximum humidity | 40 |

### OTHER FEATURES