

FEMALE TUNIC, IRIS LYOCELL , S0, HOLTEX

Ref. WMIRI07_92



09/02/2026 22:49:29 – For information purposes only, can be subject to change

FEMALE TUNIC, IRIS LYOCELL , S0, HOLTEXRef. WMIRI07_92**The eco-responsible women's medical tunic****Description**

- White women's tunic with short sleeves,
- Eco-responsible: 50% Lyocell (eucalyptus fiber) and 50% Polyester,
- Rounded neckline,
- Side opening with snap buttoning
- Nice finishes with white edging with 3 plated pockets and Kimono sleeves,
- Twill, 215g/m2,
- Available from Size 0 to Size 6.

The Holtex eco-friendly and responsible medical tunic brings extra comfort and strength to a classic poly-cotton tunic: superior quality, increased strength, easy maintenance and maximum ease. Its feminine cut is traditional and modern at the same time, perfectly corresponds to the expectations of medical and paramedical professionals: nurses, nursing assistants, doctors, nurses, physiotherapists, dental assistants, audioprosthodontists, dieticians, podiatrists, speech therapists, osteopaths, etc.

Lyocell is a natural material, made from eucalyptus wood pulp.

- Soft, lightweight, comfortable and slightly elastic,
- Easy to maintain, wrinkle resistant, it does not lint.
- Antibacterial, hypoallergenic, breathable, it absorbs moisture (up to 50% more than cotton): it is ideal for sensitive skin.
- Environmentally friendly: Eucalyptus cultivation requires much less water and pesticides than cotton.

FEMALE TUNIC, IRIS LYOCCELL , S0, HOLTEXRef. WMIRI07_92**Features**

CE Marking	NONE	Medical device	No
Bar code	3661041009420	HS Code	6211433100
Country	Tunisia		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Bleaching	Bleaching prohibited	Washing	Wash at 60°C
Model	Female model	Cleaning	Dry cleaning
Packed unit / Quantity	1	Size	Size 2