

## LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)

Ref. 519029



**CE**  
**0459**



11/02/2026 00:34:49 – For information purposes only, can be subject to change

## LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)

Ref. 519029



=====

### Description

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

#### CUFF

- Cotton CUFF, velcro fastening
- Size XXS (10-15 cm), XS (15-22 cm), S (20-28 cm, ref : 519029), M (26-34 cm, ref 519012), or L (32-42 cm), S-M-L (ref :519050)
- Available in grey, navy, red or patterns
- Supplied with carrying case

# **LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)**

**Ref. 519029**



## **Features**

|                                    |               |                           |                               |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                         | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                           | 3700446008734 | HS Code                   | 90189010                      |
| Warranty (months)                  | 36            | Country                   | France                        |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |                               |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc                            |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.663                         |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |                               |
| Minimum temperature                | -20           | Max temperature           | 70                            |
| Minimum humidity                   | 85            | Maximum humidity          | 85                            |
| <b>OTHER FEATURES</b>              |               |                           |                               |
| Packed unit / Quantity             | 1             |                           |                               |