

=====

Ref. BM103S



=====

Description

=====

 11/02/2026 00:15:43 – For information purposes only, can be subject to change

=====

Ref. BM103S

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215745	HS Code	9018901000
Warranty (months)	12		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES