

# MOBI WITH ADULT ARMBAND (M) SPENGLER, 4 COLORS TO CHOOSE FROM



Ref. 513701



**CE**  
**0459**

DISPOSITIF MEDICAL  
CLASSE  
**I<sub>m</sub>**  
MEDICAL DEVICE

GARANTIE  
**3**  
3 ANS  
GUARANTEE



📅 12/02/2026 01:15:18 – For information purposes only, can be subject to change

## MOBI WITH ADULT ARMBAND (M) SPENGLER, 4 COLORS TO CHOOSE FROM



Ref. 513701

**Elegant and colorful, the MOBI blood pressure monitor revolutionizes the aneroid blood pressure monitor thanks to its new patented connectors.**

### Description

The "LIFE LINK®" tubing allows the cuff to be changed in less than a second thanks to its magnetic connectors that ensure the connection and sealing between the pressure gauge and the tubing.

- AMBIDEXTRE : Mobi is specially designed to be ambidextrous, choose your configuration and change the position of the decompression valve, the "LIFE LINK®" connectors and the spoon.

- ERGONOMIC PEAR: The ergonomic pear is suitable for all hand morphologies and offers a volume gain of 11% on traditional pears.

- PRECISION VALVE: MOBI is equipped with a decompression valve whose shape is studied to offer the best sensitivity during decompression but also a very good ergonomics and grip.

- SHOCK RESISTANT: Due to its ABS plastic construction, its double inner shell and its elastomer rear part, MOBI is specially designed to offer great resistance to shocks and falls in everyday life.

- EASY ARMBAND with the revolutionary "LIFE LINK®" connectors, designed in TPU coating to be lighter and easy to clean.

- COLORED: Mobi exists in 4 colors: Carbon, Blueberry, Coral and Powder Pink-

Comes with an M cuff and a LIFE LINK® tubing

- Optional: 5 sizes of armbands (XS, S, M, L, XL)

# MOBI WITH ADULT ARMBAND (M) SPENGLER, 4 COLORS TO CHOOSE FROM



Ref. 513701

## Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446027964	HS Code	90189010
Warranty (months)	36	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.52
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
<b>OTHER FEATURES</b>			
Colors	Coral		

📅 12/02/2026 01:15:18 – For information purposes only, can be subject to change