

Ref. 332001



CE



16/02/2026 20:37:50 – For information purposes only, can be subject to change

---

**Ref. 332001**


## Description

Cotton cuff with highly resistant rubber bladder.  
 Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.  
 Washable at 30-40°C.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028466	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	14	Packed unit / Weight	0.31
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 16/02/2026 20:37:50 – For information purposes only, can be subject to change