

**Ref. 337846**



CE



16/02/2026 21:13:05 – For information purposes only, can be subject to change

---

**Ref. 337846**


## Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Spengler hand-held BPMs Protocol patient monitors

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446000134	HS Code	9018901000
Warranty (months)	6	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	10	Packed unit / Height (cm)	4
Packed unit / Length (cm)	13	Packed unit / Weight	0.088
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 16/02/2026 21:13:05 – For information purposes only, can be subject to change