

**Ref. 513250**



## Description

Efficiency and high-performance improved with the Clinic® version. This version of the Lian® Nano is lightweight and easy to handle and efficient thanks to its Clinic® cuff, suitable for liquid decontamination. Ideal for septic isolation or for preventing cross-contamination. Supplied in a nylon pouch.

16/02/2026 20:14:37 – For information purposes only, can be subject to change



**Ref. 513250**

## Features

CE Marking	CE0459	Medical device	No
Bar code	3700446024710	HS Code	90189010
Warranty (months)	24	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.565
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

16/02/2026 20:14:37 – For information purposes only, can be subject to change