

Ref. 513251



Description

Efficiency and high-performance improved with the Clinic® version. This version of the Lian® Nano is lightweight and easy to handle and efficient thanks to its Clinic® cuff, suitable for liquid decontamination. Ideal for septic isolation or for preventing cross-contamination. Supplied in a nylon pouch.

16/02/2026 20:41:45 – For information purposes only, can be subject to change

=====

Ref. 513251



Features

CE Marking	CE0459	Medical device	No
Bar code	3700446022693	HS Code	90189010
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.524
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

 16/02/2026 20:41:45 – For information purposes only, can be subject to change