

**SET OF 20 ADULT 23-33CM**

**Ref. DSU231065**

**SUNTECH**



CE



A single-patient adult cuff designed for reliable readings with Bayonet connector.

## Description

Set of 20 Suntech "SOFT" single-patient ambands, royal blue.  
Adult size 23-33cm, rectus bayonnette connection.

16/02/2026 20:42:07 – For information purposes only, can be subject to change

## SET OF 20 ADULT 23-33CM

Ref. DSU231065

**SUNTECH**

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014179	HS Code	9018905040
Country	China		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.319
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 16/02/2026 20:42:07 – For information purposes only, can be subject to change