

# ADU/PED/NEO NIBP TUBE (WHITE) 3M FOR NC3/NC5

COMEN

Ref. BCBRAS12

3-meter cuff extension – compatible with Comen NC3/NC5, flexible and extended reach.

## Description

This 3-meter extension connects blood pressure cuffs to Comen NC3/NC5 monitors, providing greater mobility for the patient. Ideal for clinical settings where flexibility is key.

## Features

CE Marking	NONE	Medical device	No
Bar code	3661040219332	HS Code	9018390000
Warranty (months)	6	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 17/02/2026 21:06:34 – For information purposes only, can be subject to change