

PRINTING PAPER 50MM*20M (NO LOGO) FOR NC5 (BOX OF 10 ROLLS)

Ref. BCIMPR02

High-quality ECG paper rolls for precise and continuous printing, compatible with the Comen NC5.

Description

High-quality ECG paper rolls for precise and continuous printing, compatible with the Comen NC5.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040219394	HS Code	482340000
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 17/02/2026 20:16:25 – For information purposes only, can be subject to change