

EASY II URGENCE 5 CUFFS (ADULT, XL, XXL, CHILD, INFANT, 10 À 66 CM)

Ref. DS99903



CE
0123

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
1 AN
1 YEAR
GUARANTEE



📅 17/02/2026 20:30:46 – For information purposes only, can be subject to change

EASY II URGENCE 5 CUFFS (ADULT, XL, XXL, CHILD, INFANT, 10 À 66 CM)

Ref. DS99903



Special Emergency Kit

Description

- Large screen manometer (diam 56.5cm)
 -
 - Fast and effective blood pressure
 -
 - Cuff without pocket, washable nylon with velcro closure
 -
 - Ergonomic pear with decompression valve
 -
 - Graduation from 0 to 300 mm Hg
 -
 - Shockproof ring for the dial
- Delivered with nylon storage kit and 5 armbands from 10 to 66 cm (infant, child, adult, strong arm, thigh)

NB: in accordance with the decree of 28/08/09 on the compulsory equipment of ambulances.

EASY II URGENCE 5 CUFFS (ADULT, XL, XXL, CHILD, INFANT, 10 À 66 CM)



Ref. DS99903

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3661040186924	HS Code	901890
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85

OTHER FEATURES

 17/02/2026 20:30:46 – For information purposes only, can be subject to change