

SET OF 5 SUNTECH MONOBLOC INFANT REUSABLE ARMBANDS (8-13CM) - BAYONET CONNECTOR

SUNTECH

Ref. DSU231010



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2
2 ANS
GUARANTEE

Reusable infant cuffs – size 8–13 cm, Bayonet connector, designed for gentle and accurate monitoring.

Description

Reusable nylon cuff sold by 5. Single tube, bayonet coupling.

17/02/2026 20:25:28 – For information purposes only, can be subject to change

SET OF 5 SUNTECH MONOBLOC INFANT REUSABLE ARMBANDS (8-13CM) - BAYONET CONNECTOR

SUNTECH

Ref. DSU231010

Features

CE Marking	CE	Medical device	Class I
Bar code	0120840935104138	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	9
Packed unit / Length (cm)	11	Packed unit / Weight	0.168
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 17/02/2026 20:25:28 – For information purposes only, can be subject to change