



Ref. TS50003S



CE
0459

GARANTIE
3
3 ANS
3 YEARS
GUARANTEE



📅 17/02/2026 20:26:28 – For information purposes only, can be subject to change

Ref. TS50003S



Description

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215448	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	20
Packed unit / Length (cm)	13	Packed unit / Weight	0.485
Bare Unit / Weight (kg)	0.325		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100
OTHER FEATURES			

17/02/2026 20:26:28 – For information purposes only, can be subject to change