

DISPOSABLE FLOVAC LINER 2L FOR ASPIRET, ASKIR 20, 30, C30BR

Ref. BCASPI101



18/02/2026 21:06:44 – For information purposes only, can be subject to change

DISPOSABLE FLOVAC LINER 2L FOR ASPIRET, ASKIR 20, 30, C30BR



Ref. BCASP1101

Features

CE Marking	CE0123	Medical device	No
Bar code	0000000000604	HS Code	90189084
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

18/02/2026 21:06:44 – For information purposes only, can be subject to change