

# CONTAINER FLOVAC 2L (REUSABLE) FOR ASPIRET, ASKIR 20, 30 AND C30BR



Ref. BCASP1104

## Features

CE Marking	CE0123	Medical device	No
Bar code	0000000000611	HS Code	4202121900
Country	Italy		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	14
Packed unit / Length (cm)	34	Packed unit / Weight	0.366
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

18/02/2026 21:26:14 – For information purposes only, can be subject to change